
TRU ZONE*

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Instructions / Daily Record / **ColorZone*** Tapes

⚠ Please read the instructions carefully before using.
Your **TruZone*** PFM has been designed for single patient use only.

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www.aerochamber.com

INTRODUCTION

TruZone* Peak Flow Meter (PFM) and its accompanying literature meet the most recent technical Standards of Spirometry for portable peak flow meters established by the American Thoracic Society. The **TruZone*** PFM meets the educational guidelines recommended by the National Asthma Educational Prevention Program (NAEPP) of the National Institutes of Health.

TruZone* PFM is a “one size fits all” peak flow meter, meaning that one model can be used by either children or adults. Best of all, you can personalize **TruZone*** PFM with our unique color zones as recommended by your physician within a treatment plan to meet your needs.

Just as a thermometer measures body temperature, **TruZone*** PFM measures the Peak Expiratory Flow (PEF). PEF is the fastest speed at which air is forced out from the lungs after taking a deep breath. This measurement is useful in detecting changes in airways that could signal a worsening of symptoms or an improvement in breathing function for people with respiratory conditions such as asthma or emphysema. Use the Daily Record as instructed and share this information with your physician or licensed healthcare professional at each visit.

Normal PEF numbers vary based on gender, age and height. These numbers are for reference only since each individual could have a personal best PEF number that may be higher or lower than the average. Your personal best PEF number is usually the highest PEF number that you have recorded on your Daily Record sheet over 2-3 weeks. Your personal best PEF number is your goal to reach when you use your **TruZone*** PFM. Your physician or licensed healthcare professional will advise you regarding your personal best.

Your physician or licensed healthcare professional has suggested that you use **TruZone*** PFM because it is important to measure your PEF daily as part of the overall management of your respiratory condition.

Read and follow the instructions to obtain accurate results.

The **TruZone*** PFM is intended for persons with respiratory conditions including asthma. Persons with these conditions should only use the **TruZone*** PFM and the enclosed **ColorZone*** tapes in conjunction with a treatment plan under the direction of a physician or other licensed practitioner.

If you have questions about whether **TruZone*** PFM is working correctly, or you do not fully understand the user instructions, please contact your physician.

Because your physician’s or licensed healthcare professional’s advice is important, their instructions for use must be followed carefully and consistently. They will help you understand the readings and results from your **TruZone*** PFM. The PEF is normally measured at 8 A.M. and 8 P.M., but always as recommended by a physician or licensed

healthcare professional. To insure proper technique, ask your physician or licensed healthcare professional to watch you use the **TruZone*** PFM initially.

Record your results in the Daily Record each time you use **TruZone*** PFM to identify changes, trends or patterns in your respiratory symptoms.

If you experience unusual symptoms, follow the instructions provided by your physician or licensed healthcare professional and/or contact them immediately regardless of your PEF.

Your physician or licensed healthcare professional should always be contacted when changes in your PEF indicate a worsening of your condition.

HOW TO USE YOUR **TRUZONE*** PFM

1. Prior to using your new **TruZone*** PFM make sure the yellow indicator has been reset, resting within the diamond near the mouthpiece.
2. To reset, gently shake the **TruZone*** PFM vertically until the yellow indicator is positioned directly behind the reset diamond.
3. Take as deep a breath as possible.
4. Put the **TruZone*** PFM mouthpiece into your mouth – seal your lips.
5. Blow out as hard and as fast as possible.
6. Repeat steps 2 through 5 as many times as instructed by your physician or health care professional.
7. Look at the number on the scale beside the yellow indicator.
Record this number in the space provided in your **TruZone*** PFM Daily Record.

! When you finish using your **TruZone*** PFM, be sure to store it in its original package.

HOW TO USE YOUR DAILY RECORD

The Daily Record covers a period of 35 days. Write your name on the line provided and starting with today's date, enter the dates of the next 35 days across the top of the record, in the spaces provided.

Under each date is a graph for charting PEF up to three times a day (AM, PM and one additional). This graph matches the numbers on the **TruZone*** PFM. Begin entering results in your Daily Record in the left hand column of the chart. All subsequent results should be entered sequentially to the right.

Once the yellow indicator comes to rest at your final PEF, make an "X" to record the measurement under the appropriate Date and Time heading. You have a choice of two methods when charting your peak expiratory flow results:

TruZone* PFM METHOD: Because the scales on both your **TruZone*** PFM and Daily Record are identical, it is easy to record your PEF. Simply lay your **TruZone*** PFM directly onto the Daily Record. Align the upper and lower limits of the PEF scales of both your **TruZone*** PFM and Daily Record. Now mark with an "X" below the corresponding Date and Time entry—the point where the yellow indicator rests.

Conventional Method: Simply locate the point on the Daily Record, below the Date and Time, that corresponds to the PEF you wish to record. Mark an "X" at that point.

! Always keep your Daily Record up-to-date. Additional copies of the Daily Record can be obtained by visiting www.truzone.ca

CLEANING YOUR **TRUZONE*** PFM

Your **TruZone*** PFM need only be cleaned once a week. Wash the PFM by moving gently through warm, mild, soapy water. Rinse gently and thoroughly and allow it to air dry before your next reading. Do not try to clean the inside with a brush. You can also safely clean **TruZone*** PFM using the normal cycle of your dishwasher.

CAUTION: Never disassemble your **TruZone*** PFM or place it in boiling water.

QUALITY CHECK

Always visually inspect your **TruZone*** PFM for damage. Make sure it is working properly before using. This can be done by lightly blowing into the PFM to see if the indicator moves, then resetting the indicator. A damaged or problem meter should not be used.

Note: The **TruZone*** PFM meets all accuracy and precision requirements that may be encountered under typical use. As with any pressure/flow rate monitor, the **TruZone*** PFM may be affected by extreme altitude variations or temperatures.** The meter is not significantly affected by relative humidity and should be used within a normal temperature range (20 to 35 degrees Celsius/68 to 95 degrees Fahrenheit). All measurements of the **TruZone*** PFM have been validated at Body Temperature and Pressure Saturated with water vapors (BTPS).

** Gardner, RM, et al. *Evaluation of Accuracy and Reproducibility of Peak Flow Meter at 1,400m*. Chest, 1992; 101:948-952.

STANDARDS


Your **TruZone*** PFM is designed and manufactured to meet the most recent American Thoracic Society Standards for peak flow measurements. These standards reflect the use of 26 flow-time waveforms specifically selected to evaluate peak flow meters. Accuracy is $\pm 12\%$ or 25 l/m, whichever is greater.

OPTIONAL **ColorZone*** TAPE

Our **ColorZone*** Tape system makes the use of your **TruZone*** PFM easy. These zones are based on a simple traffic light system. Your physician will instruct you on the placement and use of this unique **ColorZone*** Tape system.

ColorZone* crack-and-peel Tapes are enclosed that may help to define your individual PEF zones. Remove one piece of multicolored tape from the sheet by cutting along the dotted lines. To apply, peel the backing from the green zone. Then place the tape in the middle of the scale of the PFM so that the arrows below the word BEST line up with your best PEF as confirmed by your physician or licensed healthcare professional.

Press the green zone into place. Cut the red zone just above the reset diamond (see illustration), peel off the remaining backing and press into place.

GREEN:	80-100 percent of personal best
HIGH YELLOW:	65-80 percent of personal best
LOW YELLOW:	50-65 percent of personal best
RED:	below 50 percent of personal best

PEF may vary as your respiratory condition changes, PEF may also change as children grow. For these reasons, several **ColorZone*** Tapes are included with this packet. Used **ColorZone*** Tapes must be removed before a new one is applied. Do not replace **ColorZone*** Tapes before consulting with your physician.

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[†] Polgar G, Promadhat V. *Pulmonary Function Testing in Children: Techniques and Standards*. Philadelphia, W.B. Saunders Company, 1971.

The Lung Association acknowledges that peak flow monitoring is one useful tool to measure airway obstruction.

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Trudell Medical International*

725 Third Street, London,
Ontario, Canada, N5V 5G4
Phone: +1-519-455-7060 Fax: +1-519-455-7858
Email: customerservice@trudellmed.com
www.trudellmed.com

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Trudell Medical
Marketing Limited*

758 Third Street,
London, Ontario, Canada, N5V 5J7
Phone: +1-519-685-8800 Fax: +1-519-685-8993
Toll Free Customer Service: 1-800-265-5494
Email: customerservice@tmml.com
www.tmml.com

www.truzone.ca