

## Features & Specifications Of Medical Oxygen Cylinder - Size D



<b>Key Features</b>	
Bag	With carrying belt and compartment to store accessories
Oxygen Cylinder	Aluminium Lightweight
Oxygen Regulator	2-in-1 Pin-index Regulator and Flow meter Able to use with or without humidifier bottle Adjustable flow rate from 0 - 15 LPM

<b>Specifications</b>	
Oxygen Capacity	425 LITER
Service Pressure	2015 PSI / 139 bar
Cylinder Length	46 cm
Cylinder Outer Diameter	10 cm
Cylinder Weight (empty)	2.4 kg