

Features & Specifications Of Medical Oxygen Cylinder - E Size



Key Features	
Trolley	With rubber handle and hook to hang accessories
Oxygen Cylinder	Aluminium Lightweight
Oxygen Regulator	2-in-1 Pin-index Regulator and Flow meter Able to use with or without humidifier bottle Adjustable flow rate from 0 - 15 LPM

Specifications	
Oxygen Capacity	680 LITER
Service Pressure	2015 PSI / 139 bar
Cylinder Length	50 cm
Cylinder Outer Diameter	11 cm
Cylinder Weight (empty)	3.58 kg