

Features & Specifications Of Medical Oxygen Cylinder - Size F



Key Features	
Bag	With Rubber handle and hook to hang accessories
Oxygen Cylinder	Seamless Steel
James Regulator/Flow Meter with Humidifier connector	2-in-1 Pin-index Regulator and Flow meter Able to use with or without humidifier bottle Adjustable flow rate from 0 - 15 LPM

Specifications	
Oxygen Capacity	1360L
Service Pressure	2015 PSI / 139 bar
Cylinder Length	76 cm
Cylinder Outer Diameter	14 cm
Cylinder Weight (empty)	16 kg